

# **IDLUXHAIR**

PhotoBioModulation Treatment



# Tentech Item





# **10 LUX HAIR**

### PhotoBioModulation Treatment

10LUX-Hair is a device that increases the metabolic activity of mitochondria within cells by irradiating LEDs to the scalp, producing various effects depending on the wavelength.



# Performance





#### Various treatment modes






# Auto Lift Function

#### **10 LUX HAIR**





Moves to user's head height without using hands



# Professional Scalp Care







# Specification

Technical Specification						
Wavelength(nm) (Embedded Optical LENS)	415	633	650	830		
Max Output (mW/cm2)	30	50	20	80		
Interface	7 inch Touch screen					



Light  $\rightarrow$  Transcription factor activation  $\rightarrow$ 

#### PBM (Photo Bio Modulation)

Treatment using the principle of increasing cellular metabolic activity by light


Increased gene expression				
Protein synthesis	Increased cell activity			
Proliferation	Increased tissue repair and regeneration			
Cell Migration	Wound healing and immune response			
Anti-Inflammatory Signaling	Reduces pain and inflammation			
Anti Apoptosis	Prevent tissue damage			
Antioxidant	Prevents tissue aging			

After the initial photon absorption events, numerous signaling pathways are activated via reactive oxygen species, cyclic AMP, NO and Ca2+, leading to activation of transcription factors. These transcription factors can lead to increased expression of genes related to protein synthesis, cell migration and proliferation anti-inflammatory signaling, anti-apoptotic proteins, antioxidant enzymes.

Stem cell and progenitor cells appear to be particularly susceptible to LLLT.

\*Source: Proposed Mechanisms of Photobiomodulation or Low-Level Light Therapy. IEEE J Sel Top Quantum Electron. 2016 ; 22(3): . doi:10.1109/JSTQE.2016.2561201









#### PBM Treatment Mechanism

#### **10 LUX HAIR**

Light reacts with Cytochrome C oxydase in the mitochondrial inner membrane electron transport system to generate NO, ATP, and ROS. Therefore, it leads to an increase in gene transcription, cell regeneration, and differentiation.




2. NO(Nitric Oxide), Active oxygen: Promote cell signaling system  $\rightarrow$  Cell growth and differentiation





### Hair Treatment Mechanism





#### **10 LUX HAIR**

# Hair treatment cases using PBM



Case 1


Case 2

Case 3	3
--------	---

	Symptom	Treatment method		Result
		Target	35 people	
Case 1	Male pattern hair loss	Term	5-10 minutes, 3 times a week (6 months)	Improved hair growth in all treated areas for both men and women. Especially, It is more effective for men.
	Wavelength	655nm		
Case 2 Alopecia areata	Target	15 people	Compared to the non-irradiated area, almost 47% of patients in the laser-irradiated area indicated an increase in hair count after about 50 days.	
	Term	3 minutes (3 times a week)		
		Wavelength	600 ~1600nm	
Case 3 Male pattern hair loss		Wavelength	Blue-5J/cm2, Red-10J/cm2	Photochemical and dynamic cosmetic treatments involving laser and LED devices are effective for alopecia.
	Male pattern hair loss	Drug	Photoactive-Dr.Peel	
		Term	5 times a month	



#### **10 LUX HAIR**



#### Arndt-Schultz Curve

#### Optimization of illumination intensity and stimulation time is key.



# **ID LUX HAIR**

All information in this document is subject to copyrights and other intellectual property rights of Tentech Inc. and other licensees. No modification, reproduction or copy of any or any part of this material may be made without prior written permission from Tentech Inc. and licensees. Additional information is available upon request. The following are trademarks or registered trademarks of Tentech, Inc. in Korea and other countries.



Tentech Inc. 4F, 5F Namjeon Building, 326, Bongeunsa-ro, Gangnam-gu, Seoul, Korea TEL. 02-6954-2203 | FAX. 02-3444-4999 | www.tenlaser.com