LIPO THERA

Dual (long & normal) pulse body contouring





LIPO THERA uses non-invasive HIFU device for body slimming and fat reduction.

Mechanism of HIFU Treatment

HIFU stands for High Intensity Focused Ultrasound that non-invasively delivers high heat ($60~80~^{\circ}$ C) to the target layer (dermis layer and SMAS layer) to contract the skin layer for skin tightening and lifting.

Therefore, it causes a thermal effect only on the target layer within the skin, and does not damage the skin surface and surrounding tissue.



Cartridge

5 Transducers


A 5-dot transducer is applied, instead of the 1-dot type transducer, which makes the procedure more efficient with an irradiation area of 25cm².




1-dot type transducer




LIPO THERA



LIPO THERA Technological Features





LIPO THERA has Normal mode and Long pulse mode. By splitting one pulse into multiple pulses, Long pulse mode transfers heat more deeply, and provides continuous heat sensations with minimized pain.

(May vary from person to person)





LIPO THERA can adjust the energy level from 1 to 7, allowing customized treatment according to individual conditions.